

**TITANIUM BAR / ALL-ON-4 PROCEDURE PRESCRIPTION AND PROCEDURE GUIDE**

Required Information	
Office / Dr. Name:	Patient Name:
Office Phone:	Age: <input type="checkbox"/> M <input type="checkbox"/> F
Email:	P/U Date:
Other:	Due Date:
Preferred Contact: <input type="checkbox"/> Call <input type="checkbox"/> Email	<input type="checkbox"/> Call Me Please

Case Details			
<b>Vendor/Implant System:</b> <input type="checkbox"/> Straumann <input type="checkbox"/> Nobel <input type="checkbox"/> 3i <input type="checkbox"/> Neoss <input type="checkbox"/> Thommen <input type="checkbox"/> Other: _____	<b>Models Enclosed:</b> <input type="checkbox"/> Study Model U / L <input type="checkbox"/> Upper Model <input type="checkbox"/> Lower Model <input type="checkbox"/> Master Model <input type="checkbox"/> Upper Impression <input type="checkbox"/> Lower Impressions <input type="checkbox"/> Bite Registration	<b>Equipment Enclosed:</b> <input type="checkbox"/> Articulator <input type="checkbox"/> Articulator Box <input type="checkbox"/> Facebow <input type="checkbox"/> Meshball	<b>Parts Enclosed:</b> <input type="checkbox"/> Open Tray Impression Copings <input type="checkbox"/> 2 Sets of Analogs <input type="checkbox"/> Guide Pins (30mm) <input type="checkbox"/> Plastic Temporary Cylinders <input type="checkbox"/> Lab Screws <input type="checkbox"/> Final Screws
<b>Bar Type:</b> Fixed Removable Restoration: <input type="checkbox"/> Full Wrap Hybrid <input type="checkbox"/> Montreal Removable Restoration: <input type="checkbox"/> Attached Bar Overdenture <input type="checkbox"/> Dolder Bar <input type="checkbox"/> Paris Bar <input type="checkbox"/> Locator Attached Overdenture (No Bar)			
<b>Tooth Shade:</b> _____			

**Stages of Final Restoration:**

1	Clinical - Preliminary non-splinted impression in stock plastic tray. (Recommended 3M - ESPE Trays) LAB - Fabricate model intra-oral luting jug & custom tray	4	Clinical - Wax try-in to check bite and esthetics LAB - Scan & digitally design bar - submit for milling - reset teeth on bar & wax for final try-in
2	Clinical - Final Splinted impression in custom tray LAB - Fabricate master scanning model "Base-Bar & Wax Rim"	5	Clinical - Final try-in of milled bar (for passivity) & set up for final check of buccal & lingual contours LAB - Inject process & finish hybrid restoration
3	Clinical - VDO-records appointment LAB - Set teeth in wax for initial try-in	6	Clinical - Final inserting of prostheses!

Signature

Doctor Signature: \_\_\_\_\_ Lic#: \_\_\_\_\_